


La Marzocco Technical Bulletin #48 - NEW LINEA AV STANDARD GIGGLERS

Dear La Marzocco Family Member,

Herewith we would like to inform you about the latest modification on our LINEAS & FB/70 AV: after machine s/n 32595, ALL MACHINES are equipped with a new standard giggler, size 0,8 mm.

The giggler (# L087/GI8) will set the group flow-rate at around 700-800 ml/min, consequently allowing a better pre-infusion.



Please do not hesitate to contact me if you need further clarification.

Sincerely yours,

Enrico Wurm
After Sales Manager
La Marzocco srl