



la marzocco

handmade in florence

Preventative Maintenance Checklist

GS3 1 Year

Please fill out all sections completely

Customer Name: _____	Service Provider: _____
Phone: _____	Technician: _____
Address: _____	Other Parts Used: _____
_____	Mileage: _____
Serial Number: _____	

All Models

- | | | |
|---|---|---|
| <input type="checkbox"/> Clean and Back Flush Machine | <input type="checkbox"/> 8 mm (H.3.001) | <input type="checkbox"/> 9 mm (L105/B9) |
| <input type="checkbox"/> Replace Group Gaskets | | |
| <input type="checkbox"/> Replace Diffuser Screens (F.3.040) | | |
| <input type="checkbox"/> Replace Portafilter Baskets | Basket Type: _____ | |
| <input type="checkbox"/> Inspect Group Valve Operation | <input type="checkbox"/> Replaced <input type="checkbox"/> Acceptable | |
| <input type="checkbox"/> Test All Switches for Proper Operation | <input type="checkbox"/> Replaced <input type="checkbox"/> Acceptable | |
| <input type="checkbox"/> Test Flow Rate | Group 1: _____ | |
| <input type="checkbox"/> Test Brew Temperature Setting | Desired: _____ | Current: _____ |
| <input type="checkbox"/> Test Brew Pressure Setting | | |
| <input type="checkbox"/> Clean Auto Fill Probe | | |
| <input type="checkbox"/> Replace Steam Boiler Vacuum Breaker Valve (L193) | | |
| <input type="checkbox"/> Test Auto Fill Valve for Proper Operation | <input type="checkbox"/> Replaced <input type="checkbox"/> Acceptable | |
| <input type="checkbox"/> Replace Pressure Safety Valve (B.2.025) | | |
| <input type="checkbox"/> Inspect Drain System | | |
| <input type="checkbox"/> Test Water with La Marzocco Water Test Kit (F.9.000) | | |

If AV Model

- | | |
|--|----------------|
| <input type="checkbox"/> Test Shot Volumes | |
| <input type="checkbox"/> Test Flow Meters for Correct Resistance | Group 1: _____ |

If MP Model

- | |
|---|
| <input type="checkbox"/> Rebuild MP Valve (SA/600K) |
|---|

Notes: _____

Technician (Please Print): _____

Technician Signature: _____